



# Bee by Jenx

## Technical Specifications



	Units	Bee NM01
Age Range (approx)	mths/yrs	4mths - 5yrs
Max User Weight	kg/lbs	25 / 55
Back Angle (recline/prone)	degrees	10° - 15°
Back Pad Height (seat to top min/max)	mm/inch	270 - 360 / 10½ - 14
X-Base (footprint)	mm/inch	610 - 620 / 24 - 24½
Y-Base (footprint)	mm/inch	650 - 780 / 25½ - 30¾
Wooden Base (footprint)	mm/inch	500 - 600 / 19¾ - 23½
Chest Width (min/max)	mm/inch	150 - 270 / 6 - 10½
Hip Width (min/max)	mm/inch	130 - 250 / 5 - 9¾
Seat Height From Floor (X-Base min/max)	mm/inch	340 - 510 / 13¼ - 20
Seat Height From Floor (Y-Base min/max)	mm/inch	285 - 750 / 11¼ - 29½
Seat Height From Floor (Wooden Low Base/High Base)	mm/inch	400 - 550 / 15¾ - 21½
Seat Depth	mm/inch	110 - 350 / 4¼ - 13¾
Lower Leg Length	mm/inch	95 - 345 / 3¾ - 13½
Depth With Long Leg Conversion (min/max)	mm/inch	370 - 560 / 14½ - 22
Tilt-In-Space Angle (X and Y Base recline/prone)	degrees	30° - 10°
Tilt-In-Space Angle (Wooden Base Recline)	degrees	20°