



Supine Stander by Jenx

Technical Specifications



	Measurement	Units	Supine Stander 1	Supine Stander 2	Supine Stander 3
Age Range	approx.	mnts/yrs	9mnts - 4yrs	3yrs - 11yrs	10yrs - 18yrs
User Weight	max	kg	25	50	100
User Height	max	mm	1060	1460	1900
Chest Width	min - max	mm	130 - 310	130 - 310	220 - 440
Hip Width	min - max	mm	130 - 310	130 - 310	220 - 440
Knee Pad Width	min - max	mm	130 - 310	130 - 310	220 - 440
Tilt-in-Space	supine	degrees	10° - 90°	10° - 90°	10° - 90°
Footplate Angle	dorsiflexion - plantar flexion	degrees	0°	0°	20° - 20°
Board Height From Footplate	min - max	mm	500 - 1070	900 - 1470	1540 - 1640
Board Size	width x length	mm	480 x 940	480 x 1345	640 x 1530
Base Size	width x length	mm	575 x 660	650 x 995	750 x 1070